



DOYLE COMBO SPLINT

Description

The Doyle Combo Splint combines all the benefits of the original nasal airway splint with the expandability and comfort of sponge. This patented device is designed for use postoperatively following a septal and/or turbinate procedure. Our exclusive RHINOCELL® sponge portion expands after introduction of fluid and provides a snug yet comfortable fit. The silicone splint portion serves to stabilize the septum while providing an airway through the integral tube.

Indications for use

To maintain a nasal airway while providing septal support following surgery.

Supplied

Sterile
10 splints per dispenser carton
(1 pair per pouch, 5 pairs per carton)

Product Code

REF	
60-202	Doyle Combo Splint



FEATURES

- Dual-purpose patented design
- Superior liquid absorption and wicking
- Integral airway and suture holes

NASAL AIRWAY SPLINT

Description

The Nasal Airway Splint is designed to provide septal support and allow nasal breathing post-operatively through the integral airway. The flexible splints can be sutured through the preformed holes in the anterior tips for stabilization. Material is clear medical-grade silicone.

Indications for use

To maintain a nasal airway while providing support following surgery.

Supplied

Sterile
10 splints per dispenser carton
(1 pair per pouch, 5 pairs per carton)

Product Code

REF	
60-201	Nasal Airway Splint



FEATURES

- Clear medical grade silicone
- Integral airway
- Can be trimmed prior to insertion