



## RHINOCELL® NASAL PACKING

### Description

RHINOCELL® Nasal Packs are constructed of a unique sponge material using an exclusive formula. The small cell structure of RHINOCELL® together with a special surface treatment provides a smooth surface, which maximizes comfort while minimizing tissue ingrowth. The high tensile strength of the sponge material prevents shedding or tearing. The packs are supplied compressed for ease of insertion and will expand as fluid is introduced, giving time for accurate positioning. Available with an integral airway to facilitate nasal airflow during use. All packs are supplied with integral lanyard for securing packing.


### Indications for use

For use as a general nasal packing.

### Supplied

Sterile

10 packs per dispenser carton, individually pouched



### Nasal Packings, standard



NP-451520



NP-801520



NP-801530

### Nasal Packings, ThinPack

✓ New



TP-451020



TP-801020



TP-801030

### Nasal Packings, with airway



NP-801520A



NP-801530A

### FEATURES

- Biocompatible lint and fiber-free sponge material
- Superior liquid absorption and wicking
- Supplied compressed for ease of insertion
- Extremely soft and pliable cell structure
- Attached lanyard to secure packing

### Product Codes & Dimensions

REF	RHINOCELL® NASAL PACK	
	Dimensions (cm)	Design
NP-451520	4.5 x 1.5 x 2.0	standard
NP-801520	8.0 x 1.5 x 2.0	standard
NP-801530	8.0 x 1.5 x 3.0	standard, anatomical
TP-451020	4.5 x 1.0 x 2.0	ThinPack
TP-801020	8.0 x 1.0 x 2.0	ThinPack
TP-801030	8.0 x 1.0 x 3.0	ThinPack, anatomical
NP-801520A	8.0 x 1.5 x 2.0	with airway
NP-801530A	8.0 x 1.5 x 3.0	with airway, anatomical